


in whole brain SUVs compared to BuChE-KO controls ($p=0.037$). This significant decrease had not been observed in our previous investigations of cerebral metabolism in a BuChE expressing 5XFAD model (vs. control counterparts) at 5 months. **Conclusions:** We have previously observed variation in cerebral glucose metabolism with A β pathology in 5XFAD mice. Our current investigations suggest that BuChE may also be a modulator of cerebral glucose metabolism in AD, whereby the presence/absence of BuChE could significantly impact not only pathology but also brain function over the course of AD. Comparison of these 5XFAD/BuChE-KO brain metabolism results with those in BuChE expressing 5XFAD counterparts is currently underway. Furthermore, regional assessment of brain metabolism and corroborating neuropathology in these animals may further implicate BuChE in the progression of AD and may highlight the potential role of targeted BuChE imaging approaches for AD diagnostics.

IC-P-023 **PILOT STUDY OF ALTERED COPPER METABOLISM AS A BIOMARKER FOR EARLY DIAGNOSIS OF ALZHEIMER'S DISEASE WITH $^{64}\text{CuCl}_2$ -PET/CT**

Fangyu Peng¹, Otto Muzik², ¹University of Texas Southwestern Medical Center, Dallas, TX, USA; ²Wayne State University, Detroit, MI, USA. Contact e-mail: fangyu.peng@utsouthwestern.edu

Background: Copper is a trace element required for development and normal function of human brains. Emerging body of evidence suggests the role of copper in pathogenesis of Alzheimer's disease. The aim of this study was to explore the potential of altered copper metabolism as a biomarker for early diagnosis of Alzheimer's disease with positron emission tomography/computed tomography (PET/CT) using copper-64 chloride ($^{64}\text{CuCl}_2$) as a radioactive tracer ($^{64}\text{CuCl}_2$ -PET/CT). **Methods:** APPSWE transgenic mice (N=4, 14 weeks old.), a mouse model of Alzheimer's disease, were subjected to PET/CT after intravenous injection of copper-64 chloride ($^{64}\text{CuCl}_2$) as a tracer, using a small animal PET/CT scanner. A group of wild type C57BL/6 mice (N=4, 13 weeks old) and another group of *Atp7b*^{-/-} knockout mice (N=4, 6 to 7 weeks old), a mouse model of Wilson's disease, were used as controls. PET quantitative analysis was performed to compare ^{64}Cu



uptake in the brains of the APPSWE transgenic mice with the ^{64}Cu uptake in the brains of C57BL/6 and *Atp7b*^{-/-} KO mice, respectively. **Results:** Increased ^{64}Cu uptake was detected in the brains of the APPSWE transgenic mice, compared with the ^{64}Cu uptake in the brains of the C57BL/6 mice and the *Atp7b*^{-/-} KO mice, respectively. In addition to increased ^{64}Cu uptake in the cortex, large increase of ^{64}Cu uptake was also detected in the regions of olfactory bulb and caudate of the APPSWE transgenic mice. Furthermore, cerebral ^{64}Cu uptake in the brains of the *Atp7b*^{-/-} KO mice was found to be lower than the cerebral ^{64}Cu uptake in both the C57BL/6 mice and the APPSWE transgenic mice. Decrease of cerebral ^{64}Cu uptake in the *Atp7b*^{-/-} KO mice was likely secondary to metabolic trapping of ^{64}Cu in the liver of *Atp7b*^{-/-} KO mice as visualized on the PET/CT images. **Conclusions:** Increased ^{64}Cu uptake was detected in the brains of APPSWE transgenic mice, compared with the ^{64}Cu uptake in the brains of C57BL/6 mice and *Atp7b*^{-/-} KO mice, respectively. The findings support further investigation of altered copper metabolism as a biomarker for early diagnosis of Alzheimer's disease with PET/CT using $^{64}\text{CuCl}_2$ as a radioactive tracer ($^{64}\text{CuCl}_2$ -PET/CT).

IC-P-024 **A NOVEL POSITRON EMISSION TOMOGRAPHY CONTRAST AGENT TARGETING CATHEPSIN D SHOWS PREFERENTIAL *IN VIVO* RETENTION IN AN ALZHEIMER'S DISEASE MOUSE MODEL**

Jonatan A. Snir^{1,2}, Mojmir Suchy^{1,2}, Geron A. Bindseil², Blaine A. Chronik², Robert H.E. Hudson², Stephen H. Pasternak^{1,2}, Robert Bartha^{1,2}, ¹Robarts Research Institute, London, ON, Canada; ²Western University, London, ON, Canada. Contact e-mail: jnsnir@robarts.ca

Background: Early detection of Alzheimer's disease (AD) pathology remains a serious challenge for both diagnosis and development of treatment. Cathepsin D (CatD), a lysosomal aspartyl protease, is over-expressed in AD and therefore is a potential biomarker. Previously, we introduced a novel Contrast Agent (CA) that was preferentially taken up by CatD over-expressing cells (*in-vitro*) and able to transverse the BBB in mice (*ex-vivo*). We have also found that a Near-Infrared-labeled version of this CA demonstrates prolonged *in-vivo* retention in the brain of a transgenic (Tg) mouse model of AD at 12 months compared to age matched wild type controls. Here, we present the performance of a CA labeled with ^{68}Ga evaluated by micro Positron Emission Tomography (microPET). **Methods:** The CA consists of a Cell Penetrating Peptide (CPP; the Tat peptide from HIV-1), attached to a CatD cleavage sequence followed by a ^{68}Ga labeled DOTA chelator flanked by a fluorescent dye. The CPP allows the agent to cross the blood brain barrier bidirectionally. In the presence of elevated levels of CatD, cleavage of the CatD site removes the CPP resulting in prolonged retention of the imaging moiety. For this study, Tg AD mice (N=8, 5XFAD model) and non-Tg age matched littermates (N=8) at 2 and 4.5 months of age received an intravenous tail vein CA injection of ~12 MBq of CA under isoflurane anesthesia. Mice were scanned for 3 hours using the Inveon preclinical microPET system (Siemens Medical Solutions, Knoxville TN, USA). Regions of interest were identified in reconstructed images and were used to measure the uptake and washout of the CA in the brain, liver, kidneys and bladder. **Results:** The Tg mice demonstrated significantly greater uptake ($p<0.05$) of the CA in the brain in the first two hours following injection at 4.5 months but not 2 months of age compared to controls.